



Plaur Antibody, Biotin conjugated

Product Code	CSB-PA018122LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49616
Immunogen	Recombinant Rat Urokinase plasminogen activator surface receptor protein (25-299AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Urokinase plasminogen activator surface receptor (U-PAR) (uPAR) (CD_antigen: CD87), Plaur
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Plaur
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.