



TSHR Antibody, Biotin conjugated

Product Code	CSB-PA11189D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16473
Immunogen	Recombinant Human Thyrotropin receptor protein (21-413AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyrotropin receptor (Thyroid-stimulating hormone receptor) (TSH-R), TSHR, LGR3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TSHR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.