



ptxA Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA356423LB01Bua |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P04977 |
| Immunogen | Recombinant Bordetella pertussis Pertussis toxin subunit 1 protein (35-269AA) |
| Raised In | Rabbit |
| Species Reactivity | Bordetella pertussis |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Pertussis toxin subunit 1 (PTX S1) (Islet-activating protein S1) (IAP S1) (NAD-dependent ADP-ribosyltransferase) (EC 2.4.2.-), ptxA |
| Immunogen Species | Bordetella pertussis |
| Research Area | Others |
| Target Names | ptxA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |