




ptxA Antibody

Product Code	CSB-PA356423LA01Bua
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04977
Immunogen	Recombinant Bordetella pertussis Pertussis toxin subunit 1 protein (35-269AA)
Raised In	Rabbit
Species Reactivity	Bordetella pertussis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pertussis toxin subunit 1 (PTX S1) (Islet-activating protein S1) (IAP S1) (NAD-dependent ADP-ribosyltransferase) (EC 2.4.2.-), ptxA
Immunogen Species	Bordetella pertussis
Research Area	Others
Target Names	ptxA

Image



Western blot

All lanes: Bordella pertussis pertussis toxin subunit 1 antibody at 2µg/ml + recombinant Bordella pertussis pertussis toxin subunit 1 100ng

Secondary

Goat polyclonal to rabbit IgG at 1/1000 dilution
Predicted band size: 36 kDa
Observed band size: 36 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.