



# SOSTDC1 Antibody, FITC conjugated

|                            |                                                                                                                                                           |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA022416EC01HU                                                                                                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                             |
| <b>Uniprot No.</b>         | Q6X4U4                                                                                                                                                    |
| <b>Immunogen</b>           | Recombinant Human Sclerostin domain-containing protein 1 protein (24-206AA)                                                                               |
| <b>Raised In</b>           | Rabbit                                                                                                                                                    |
| <b>Species Reactivity</b>  | Human                                                                                                                                                     |
| <b>Form</b>                | Liquid                                                                                                                                                    |
| <b>Conjugate</b>           | FITC                                                                                                                                                      |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                                                          |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                                  |
| <b>Isotype</b>             | IgG                                                                                                                                                       |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                |
| <b>Alias</b>               | Sclerostin domain-containing protein 1 (Ectodermal BMP inhibitor) (Ectodin)<br>(Uterine sensitization-associated gene 1 protein) (USAG-1), SOSTDC1, USAG1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                      |
| <b>Research Area</b>       | Signal Transduction                                                                                                                                       |
| <b>Target Names</b>        | SOSTDC1                                                                                                                                                   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                               |