



SOSTDC1 Antibody

Product Code	CSB-PA022416EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6X4U4
Immunogen	Recombinant Human Sclerostin domain-containing protein 1 protein (24-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sclerostin domain-containing protein 1 (Ectodermal BMP inhibitor) (Ectodin) (Uterine sensitization-associated gene 1 protein) (USAG-1), SOSTDC1, USAG1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SOSTDC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.