




NAPRT Antibody

Product Code	CSB-PA015451EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XQN6
Immunogen	Recombinant Human Nicotinate phosphoribosyltransferase protein (1-310AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NAPRT

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, Mouse liver tissue, Mouse kidney tissue
 All lanes: NAPRT antibody at 4µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 58, 61, 57 kDa
 Observed band size: 58 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.