



MRJP1 Antibody, HRP conjugated

Product Code	CSB-PA522725EB01DNK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O18330
Immunogen	Recombinant Apis mellifera Major royal jelly protein 1 protein (20-432AA)
Raised In	Rabbit
Species Reactivity	Apis mellifera
Tested Applications	ELISA
Form	liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major royal jelly protein 1 (MRJP-1) (56-kDa protein 4) (p56kP-4) (Bee-milk protein) (Royalactin) [Cleaved into: Jellein-1 (Jelleine-I); Jellein-2 (Jelleine-II); Jellein-4 (Jelleine-IV)], MRJP1
Immunogen Species	Apis mellifera (Honeybee)
Research Area	Others
Target Names	MRJP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.