




Ca1 Antibody

Product Code	CSB-PA004364HA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B0BNN3
Immunogen	Recombinant Rat Carbonic anhydrase 1 protein (2-261AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carbonic anhydrase 1 (EC 4.2.1.1) (Carbonate dehydratase I) (Carbonic anhydrase I) (CA-I), Ca1, Car1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Ca1

Image



Western blot
 All lanes: Ca1 antibody at 2µg/ml + Rat spleen tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 29 kDa
 Observed band size: 29 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.