



H2-K1 Antibody

Product Code	CSB-PA355780EA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01902
Immunogen	Recombinant Mouse H-2 class I histocompatibility antigen, K-D alpha chain protein (22-305AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	H-2 class I histocompatibility antigen, K-D alpha chain (H-2K(D)), H2-K1, H2-K
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	H2-K1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.