



GALNS Antibody, Biotin conjugated

Product Code	CSB-PA009202ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34059
Immunogen	Recombinant Human N-acetylgalactosamine-6-sulfatase protein (27-522AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acetylgalactosamine-6-sulfatase (EC 3.1.6.4) (Chondroitinsulfatase) (Chondroitinase) (Galactose-6-sulfate sulfatase) (GalN6S) (N-acetylgalactosamine-6-sulfate sulfatase) (GalNAc6S sulfatase), GALNS
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GALNS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.