



IFNA2 Antibody, Biotin conjugated

Product Code	CSB-PA06504D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01563
Immunogen	Recombinant Human Interferon alpha-2 protein (24-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon alpha-2 (IFN-alpha-2) (Interferon alpha-A) (LeIF A), IFNA2, IFNA2A IFNA2B IFNA2C
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFNA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.