



IDO1 Antibody, Biotin conjugated

Product Code	CSB-PA010996HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14902
Immunogen	Recombinant Human Indoleamine 2,3-dioxygenase 1 protein (1-403AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Indoleamine 2,3-dioxygenase 1 (IDO-1) (EC 1.13.11.52) (Indoleamine-pyrrole 2,3-dioxygenase), IDO1, IDO INDO
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IDO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.