



PARN Antibody, HRP conjugated

Product Code	CSB-PA017456LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95453
Immunogen	Recombinant Human Poly(A)-specific ribonuclease PARN protein (1-639AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Poly(A)-specific ribonuclease PARN (EC 3.1.13.4) (Deadenylating nuclease) (Deadenylation nuclease) (Polyadenylate-specific ribonuclease), PARN, DAN
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PARN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.