



XPC Antibody

Product Code	CSB-PA026217LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01831
Immunogen	Recombinant Human DNA repair protein complementing XP-C cells protein (496-734AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA repair protein complementing XP-C cells (Xeroderma pigmentosum group C-complementing protein) (p125), XPC, XPCC
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	XPC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.