



FGFR3 Antibody, FITC conjugated

Product Code	CSB-PA008646LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22607
Immunogen	Recombinant Human Fibroblast growth factor receptor 3 protein (23-375AA)
Raised In	Rabbit
Species Reactivity	Human
Form	liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor receptor 3 (FGFR-3) (EC 2.7.10.1) (CD antigen CD333), FGFR3, JTK4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FGFR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.