




HSPA6 Antibody

Product Code	CSB-PA863082HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17066
Immunogen	Recombinant Human Heat shock 70 kDa protein 6 protein (1-637AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock 70 kDa protein 6 (Heat shock 70 kDa protein B'), HSPA6, HSP70B'
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPA6

Image



Western blot
 All lanes: HSPA6 antibody at 2µg/ml
 Lane 1: MCF-7 whole cell lysate
 Lane 2: A431 whole cell lysate
 Lane 3: Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 72 kDa
 Observed band size: 72 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.