



FCGRT Antibody, FITC conjugated

Product Code	CSB-PA008545EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55899
Immunogen	Recombinant Human IgG receptor FcRn large subunit p51 protein (24-297AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	IgG receptor FcRn large subunit p51 (FcRn) (IgG Fc fragment receptor transporter alpha chain) (Neonatal Fc receptor), FCGRT, FCRN
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCGRT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.