




Mmp7 Antibody

Product Code	CSB-PA014677EA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50280
Immunogen	Recombinant Rat Matrilysin protein (98-267AA)
Raised In	Rabbit
Species Reactivity	Rat, Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrilysin (EC 3.4.24.23) (Matrin) (Matrix metalloproteinase-7) (MMP-7) (Pump-1 protease) (Uterine metalloproteinase), Mmp7, Mmp-7
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Mmp7

Image



Western blot

All lanes: Mmp7 antibody at 2µg/ml

Lane 1: Colo320 whole cell lysate

Lane 2: NIH/3T3 whole cell lysate

Lane 3: Mouse lung tissue

Lane 4: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.