



CEACAM4 Antibody

Product Code	CSB-PA005164LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75871
Immunogen	Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 4 protein (36-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carcinoembryonic antigen-related cell adhesion molecule 4 (Carcinoembryonic antigen CGM7) (Non-specific cross-reacting antigen W236), CEACAM4, CGM7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CEACAM4