



Cd163 Antibody, Biotin conjugated

Product Code	CSB-PA650452LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2VLH6
Immunogen	Recombinant Mouse Scavenger receptor cysteine-rich type 1 protein M130 protein (86-365AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Scavenger receptor cysteine-rich type 1 protein M130 (CD antigen CD163) [Cleaved into: Soluble CD163 (sCD163)], Cd163, M130
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cd163