



lapB Antibody, FITC conjugated

Product Code	CSB-PA359577HC01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AB58
Immunogen	Recombinant Escherichia coli Lipopolysaccharide assembly protein B protein (17-389AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipopolysaccharide assembly protein B (Lipopolysaccharide regulatory protein), lapB, yciM
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	lapB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.