




lapB Antibody

Product Code	CSB-PA359577HA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AB58
Immunogen	Recombinant Escherichia coli Lipopolysaccharide assembly protein B protein (17-389AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	lapB

Image



Western blot

All lanes: lapB antibody at 1: 1000

Lane1: DH5α whole cell lysate at 10μg

Lane2: DH5α whole cell lysate at 1μg

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 45 kDa

Observed band size: 45 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.