



glpD Antibody, HRP conjugated

Product Code	CSB-PA320174HB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13035
Immunogen	Recombinant Escherichia coli Aerobic glycerol-3-phosphate dehydrogenase protein (1-501AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aerobic glycerol-3-phosphate dehydrogenase (EC 1.1.5.3), glpD, glyD
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	glpD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.