




RBP1 Antibody

Product Code	CSB-PA019480YA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09455
Immunogen	Recombinant Human Retinol-binding protein 1 protein (70-112AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinol-binding protein 1 (Cellular retinol-binding protein) (CRBP) (Cellular retinol-binding protein I) (CRBP-I), RBP1, CRBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RBP1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA019480YA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.