



LGB Antibody, HRP conjugated

Product Code	CSB-PA00280F0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02754
Immunogen	Bovine milk Beta-lactoglobulin (Native Protein)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-lactoglobulin (Beta-LG) (allergen Bos d 5), LGB
Immunogen Species	Bos taurus (Bovine)
Research Area	Cell Biology
Target Names	LGB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.