




α 1-acidglycoprotein Antibody

Product Code	CSB-PA00330E0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Human α 1-acidglycoprotein (Native Protein)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	α 1-AGP
Immunogen Species	Homo sapiens (Human)
Research Area	Others

Image



Western blot

All lanes: α 1-acidglycoprotein antibody at 2 μ g/ml
+ Human serum at 1: 50
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 45 kDa
Observed band size: 45 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.