



GPI Antibody, Biotin conjugated

Product Code	CSB-PA00367D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06744
Immunogen	Recombinant Human Glucose-6-phosphate isomerase protein (2-554AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucose-6-phosphate isomerase (GPI) (EC 5.3.1.9) (Autocrine motility factor) (AMF) (Neuroleukin) (NLK) (Phosphoglucose isomerase) (PGI) (Phosphohexose isomerase) (PHI) (Sperm antigen 36) (SA-36), GPI
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPI
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.