



Wheat Gliadin Antibody, HRP conjugated

Product Code	CSB-PA00450F0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Wheat Gliadin (Native Protein)
Raised In	Rabbit
Species Reactivity	Plant
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Wheat Gliadin
Immunogen Species	Plant
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.