





RHOA Antibody

Product Code	CSB-PA00875A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61586
Immunogen	Recombinant Human Transforming protein RhoA protein (1-191AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming protein RhoA (Rho cDNA clone 12) (h12), RHOA, ARH12 ARHA RHO12
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RHOA

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA00875A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA00875A0Rb at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.