



MYH9 Antibody, HRP conjugated

Product Code	CSB-PA01405B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35579
Immunogen	Recombinant Human Myosin-9 protein (2-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-9 (Cellular myosin heavy chain, type A) (Myosin heavy chain 9) (Myosin heavy chain, non-muscle IIa) (Non-muscle myosin heavy chain A) (NMMHC-A) (Non-muscle myosin heavy chain IIa) (NMMHC II-a) (NMMHC-IIA), MYH9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYH9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.