





GATD3A Antibody

Product Code	CSB-PA01814A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DPI2
Immunogen	Recombinant Human Glutamine amidotransferase-like class 1 domain-containing protein 3A, mitochondrial protein (43-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamine amidotransferase-like class 1 domain-containing protein 3A, mitochondrial, GATD3A, C21orf33
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GATD3A


Image



Western blot
 All lanes: GATD3A antibody at 2μg/ml
 Lane 1: EC109 whole cell lysate
 Lane 2: 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/15000 dilution
 Predicted band size: 29, 25 kDa
 Observed band size: 29 kDa



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA01814A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA01814A0Rb at dilution of 1:100