



FTH1 Antibody, Biotin conjugated

Product Code	CSB-PA02089D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02794
Immunogen	Recombinant Human Ferritin heavy chain protein (2-183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ferritin heavy chain (Ferritin H subunit) (EC 1.16.3.1) (Cell proliferation-inducing gene 15 protein) [Cleaved into: Ferritin heavy chain, N-terminally processed], FTH1, FTH FTHL6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FTH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.