



ZBED1 Antibody, Biotin conjugated

Product Code	CSB-PA02715D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O96006
Immunogen	Recombinant Human Zinc finger BED domain-containing protein 1 protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger BED domain-containing protein 1 (Putative Ac-like transposable element) (dREF homolog), ZBED1, ALTE DREF KIAA0785 TRAMP
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZBED1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.