



# ZNF384 Antibody, HRP conjugated

Product Code	CSB-PA03455B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TF68
Immunogen	Recombinant Human Zinc finger protein 384 protein (1-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 384 (CAG repeat protein 1) (CAS-interacting zinc finger protein) (Nuclear matrix transcription factor 4) (Nuclear matrix protein 4) (Trinucleotide repeat-containing gene 1 protein), ZNF384, CAGH1 CIZ NMP4 TNRC1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF384