



ZNF384 Antibody, Biotin conjugated

Product Code	CSB-PA03455D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TF68
Immunogen	Recombinant Human Zinc finger protein 384 protein (1-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 384 (CAG repeat protein 1) (CAS-interacting zinc finger protein) (Nuclear matrix transcription factor 4) (Nuclear matrix protein 4) (Trinucleotide repeat-containing gene 1 protein), ZNF384, CAGH1 CIZ NMP4 TNRC1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF384
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.