



RAD23A Antibody, FITC conjugated

Product Code	CSB-PA03774C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54725
Immunogen	Recombinant Human UV excision repair protein RAD23 homolog A protein (2-201AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UV excision repair protein RAD23 homolog A (HR23A) (hHR23A), RAD23A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RAD23A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.