



virD2 Antibody, Biotin conjugated

Product Code	CSB-PA03944D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18592
Immunogen	Recombinant Agrobacterium tumefaciens T-DNA border endonuclease VirD2 protein (1-447AA)
Raised In	Rabbit
Species Reactivity	Agrobacterium tumefaciens
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-DNA border endonuclease VirD2 (EC 3.1.-.-), virD2
Immunogen Species	Agrobacterium tumefaciens
Research Area	Others
Target Names	virD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.