



GPR161 Antibody, FITC conjugated

Product Code	CSB-PA04114C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6U8
Immunogen	Recombinant Human G-protein coupled receptor 161 protein (1-28AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein coupled receptor 161 (G-protein coupled receptor RE2), GPR161
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GPR161
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.