




CNPY2 Antibody

Product Code	CSB-PA05034A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2B0
Immunogen	Recombinant Human Protein canopy homolog 2 protein (22-182AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein canopy homolog 2 (MIR-interacting saposin-like protein) (Putative secreted protein Zsig9) (Transmembrane protein 4), CNPY2, MSAP TMEM4 ZSIG9
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CNPY2

Image



Western blot


All lanes: CNPY2 antibody at 2µg/ml + 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 21, 10 kDa

Observed band size: 47 kDa



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA05034A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.