




ANAPC5 Antibody

Product Code	CSB-PA05104A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJX4
Immunogen	Recombinant Human Anaphase-promoting complex subunit 5 protein (2-232AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anaphase-promoting complex subunit 5 (APC5) (Cyclosome subunit 5), ANAPC5, APC5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ANAPC5

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA05104A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.