





CTDNEP1 Antibody

Product Code	CSB-PA05415A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95476
Immunogen	Recombinant Human CTD nuclear envelope phosphatase 1 protein (3-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CTD nuclear envelope phosphatase 1 (EC 3.1.3.16) (Serine/threonine-protein phosphatase dullard), CTDNEP1, DULLARD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CTDNEP1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA05415A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA05415A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.