



IL12B Antibody, HRP conjugated

Product Code	CSB-PA06615B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29460
Immunogen	Recombinant Human Interleukin-12 subunit beta protein (23-328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-12 subunit beta (IL-12B) (Cytotoxic lymphocyte maturation factor 40 kDa subunit) (CLMF p40) (IL-12 subunit p40) (NK cell stimulatory factor chain 2) (NKSF2), IL12B, NKSF2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL12B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.