



IL6 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA06757C0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P05231 |
| Immunogen | Recombinant Human Interleukin-6 protein (30-212AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Interleukin-6 (IL-6) (B-cell stimulatory factor 2) (BSF-2) (CTL differentiation factor) (CDF) (Hybridoma growth factor) (Interferon beta-2) (IFN-beta-2), IL6, IFNB2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | IL6 |