



LCN2 Antibody, Biotin conjugated

Product Code	CSB-PA06837D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80188
Immunogen	Recombinant Human Neutrophil gelatinase-associated lipocalin protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutrophil gelatinase-associated lipocalin (NGAL) (25 kDa alpha-2-microglobulin-related subunit of MMP-9) (Lipocalin-2) (Oncogene 24p3) (Siderocalin LCN2) (p25), LCN2, HNL NGAL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LCN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.