



LCN2 Antibody, Biotin conjugated

Product Code	CSB-PA06837D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80188
Immunogen	Recombinant Human Neutrophil gelatinase-associated lipocalin protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutrophil gelatinase-associated lipocalin (NGAL) (25 kDa alpha-2-microglobulin-related subunit of MMP-9) (Lipocalin-2) (Oncogene 24p3) (Siderocalin LCN2) (p25), LCN2, HNL NGAL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LCN2