



SPP1 Antibody, Biotin conjugated

Product Code	CSB-PA06904D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10451
Immunogen	Recombinant Human Osteopontin protein (17-302AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Osteopontin (Bone sialoprotein 1) (Nephropontin) (Secreted phosphoprotein 1) (SPP-1) (Urinary stone protein) (Uropontin), SPP1, BNSP OPN
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SPP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.