



TGFB1 Antibody, Biotin conjugated

Product Code	CSB-PA07309D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01137
Immunogen	Recombinant Human Transforming growth factor beta-1 proprotein protein (281-390AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-1 (TGF-beta-1) [Cleaved into: Latency-associated peptide (LAP)], TGFB1, TGFB
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TGFB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.