




TNFRSF10D Antibody

Product Code	CSB-PA08209A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBN6
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 10D protein (58-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 10D (Decoy receptor 2) (DcR2) (TNF-related apoptosis-inducing ligand receptor 4) (TRAIL receptor 4) (TRAIL-R4) (TRAIL receptor with a truncated death domain) (CD antigen CD264), TNFRSF10D, DCR2 TRAILR4 TRUNDD
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TNFRSF10D

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA08209A0Rb at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.