



Cxcl5 Antibody, HRP conjugated

Product Code	CSB-PA09257B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50228
Immunogen	Recombinant Mouse C-X-C motif chemokine 5 protein (46-117AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 5 (Cytokine LIX) (Small-inducible cytokine B5) [Cleaved into: GCP-2(1-78); GCP-2(9-78)], Cxcl5, Scyb5
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cxcl5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.